

Anti-Ereg Rabbit Polyclonal Antibody For Wb, Elisa - O14944

Item Number: CM0024051



Introduction

Rabbit polyclonal antibody targeting human EREG (proepiregulin) for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for EGFR signaling, cancer, and inflammation research applications.

[Learn More](#)

Parameter	Value
Product Name	EREG Rabbit pAb
Remarks/Alias	ER, Ep, EPR, EREG
Species	Human
Gene ID (Human)	2069
Gene ID	2069
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-169 of human EREG (NP_001423.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14944
Protein Weight	19kDa
Shipping	Ice bag