


Anti-Erf3A/Gspt1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15170

Item Number: CM0012525



Introduction

Rabbit monoclonal antibody against eRF3A/GSPT1, reactive in Human and Rat, for WB, IF/ICC, IF-P, IHC-P, ELISA. Essential for translation termination and nonsense-mediated mRNA decay.

[Learn More](#)

Parameter	Value
Product Name	eRF3A/GSPT1 Rabbit mAb
Remarks/Alias	GST1, ETF3A, eRF3a, 551G9.2
Species	Human
GeneID (Human)	2935
GeneID	2935
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of eRF3A/GSPT1 (NP_001123479.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P15170
Protein Weight	56kDa
Shipping	Ice bag