


Anti-Erg Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11308

Item Number: CM0007666



Introduction

High-performance rabbit monoclonal antibody targeting human ERG (SWISS P11308). Validated for Western blot, IHC-P, and ELISA. Ideal for cancer research, transcription regulation studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ERG Rabbit mAb
Remarks/Alias	p55; erg-3; ERG
Species	Human
GeneID (Human)	2078
GeneID	2078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 379-479 of human ERG. (NP_891548.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11308
Protein Weight	54 kDa
Shipping	Ice bag