


Anti-Erg Rabbit Polyclonal Antibody For Wb, Elisa - P11308

Item Number: CM0015110



Introduction

High-quality rabbit polyclonal antibody against human ERG, validated for Western blot and ELISA. Suitable for transcriptional regulation and cancer research studies. - P11308

[Learn More](#)

Parameter	Value
Product Name	ERG Rabbit pAb
Remarks/Alias	p55; erg-3; ERG
Species	Human
Gene ID (Human)	2078
Gene ID	2078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ERG (NP_891548.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P11308
Protein Weight	54kDa
Shipping	Ice bag