

Anti-Ergic2 Polyclonal Antibody For Wb - Q96Rq1

Item Number: CM0027843

Introduction

Rabbit polyclonal antibody against human ERGIC2, validated for WB and ELISA. Cross-reacts with mouse. Targets ER-Golgi transport protein, 43 kDa. Ideal for intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	ERGIC2 Rabbit pAb
Remarks/Alias	PTX1; CDA14; Erv41; cd002; ERGIC2
Species	Human
GeneID (Human)	51290
Gene ID	51290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-310 of human ERGIC2 (NP_057654.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RQ1
Protein Weight	43kDa
Shipping	Ice bag