

Anti-Erk1/2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28482/P27361

Item Number: CM0010941



Introduction

CamelBio's rabbit monoclonal antibody against ERK1/2 (MAPK1/3). Validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying MAPK/ERK signaling in cell growth and differentiation.

[Learn More](#)

Parameter	Value
Product Name	ERK1/2 Rabbit mAb
Remarks/Alias	ERK; ERK-2; ERK2; ERT1; MAPK2; P42MAPK; PRKM1; PRKM2; p38; p40; p41; p41mapk; p42-MAPK; 5594/5595; ERK1/2
Species	Human
GeneID (Human)	5594/5595
GeneID	5594/5595
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-360 of human ERK2 (P28482).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	42 kDa/44 kDa
Shipping	Ice bag