

Anti-Erk1/Erk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P28482 / P27361

Item Number: CM0022556



Introduction

High-quality rabbit polyclonal antibody against human ERK1/ERK2 (MAPK3/MAPK1). Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying MAPK/ERK signaling in cell growth, differentiation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ERK1 / ERK2 Rabbit pAb
Remarks/Alias	MAPK1/MAPK3; ERK1 / ERK2
Species	Human
GeneID (Human)	5594/5595
GeneID	5594/5595
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	36kDa/41kDa/38kDa/40kDa/43kDa
Shipping	Ice bag