


Anti-Erk2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P28482

Item Number: CM0003267



Introduction

Rabbit monoclonal antibody targeting ERK2 (MAPK1). Validated for WB, IHC-P, ELISA with cross-reactivity in human, mouse, rat. Ideal for studying the MAPK/ERK signaling pathway in proliferation, differentiation, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	ERK2 Rabbit mAb
Remarks/Alias	ERK; p38; p40; p41; ERK2; ERT1; NS13; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK
Species	Human
GeneID (Human)	5594
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ERK2/ERK1 (P27361/P28482)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482
Protein Weight	41 kDa
Shipping	Ice bag