

Anti-Erk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28482

Item Number: CM0014390



Introduction

KO-validated rabbit polyclonal antibody targeting human ERK2 (MAPK1) for use in WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. An essential tool for MAPK/ERK pathway and cell proliferation investigations.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ERK2 Rabbit pAb
Remarks/Alias	ERK; p38; p40; p41; ERK2; ERT1; NS13; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK; K2
Species	Human
GeneID (Human)	5594
GeneID	5594
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human ERK2 (NP_002736.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482
Protein Weight	41kDa
Shipping	Ice bag