

Anti-Erk5 Rabbit Monoclonal Antibody For Wb, Elisa - Q13164

Item Number: CM0011501



Introduction

ERK5 Rabbit Monoclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for studying MAPK signaling, cell proliferation, differentiation, apoptosis, and cardiac biology.

[Learn More](#)

Parameter	Value
Product Name	ERK5 Rabbit mAb
Remarks/Alias	BMK1, ERK4, ERK5, PRKM7
Species	Human
GeneID	5598
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13164
Protein Weight	88kDa
Shipping	Ice bag