

Anti-Erk5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13164

Item Number: CM0018899



Introduction

Rabbit polyclonal antibody against human ERK5 (MAPK7/BMK1), validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Targets MAPK signaling pathway, cell proliferation, and differentiation.

[Learn More](#)

Parameter	Value
Product Name	ERK5 Rabbit pAb
Remarks/Alias	BMK1; ERK4; ERK5; PRKM7
Species	Human
Gene ID (Human)	5598
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13164
Protein Weight	88kDa
Shipping	Ice bag