


Anti-Erlin-2 Monoclonal Antibody For Wb, Ihc-P, Elisa - O94905

Item Number: CM0002545



Introduction

Rabbit monoclonal antibody against human Erlin-2 (ERLIN2, SPFH2). Validated for Western blot, IHC-P, and ELISA. Detects human, mouse, rat. Ideal for ERAD and cholesterol metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Erlin-2 Rabbit mAb
Remarks/Alias	NET32; SPFH2; SPG18; C8orf2; Erlin-2
Species	Human
GeneID (Human)	11160
GeneID	11160
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Erlin-2 (O94905).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94905
Protein Weight	38kDa
Shipping	Ice bag