


Anti-Erlin1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - O75477

Item Number: CM0013055



Introduction

Rabbit polyclonal antibody against human ERLIN1, suitable for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Designed to support research on ER-associated degradation and cholesterol metabolism.

[Learn More](#)

Parameter	Value
Product Name	ERLIN1 Rabbit pAb
Remarks/Alias	KE04, KEO4, SPFH1, SPG62, Erlin-1, C10orf69, ERLIN1
Species	Human
GeneID	10613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 259-348 of human ERLIN1 (NP_006450.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75477
Protein Weight	39kDa
Shipping	Ice bag