

Anti-Ermap Rabbit Polyclonal Antibody For Wb, Elisa - Q96PI5

Item Number: CM0015473

Introduction

Anti-ERMAP rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human and mouse ERMAP, a 53kDa erythroid membrane protein involved in cell adhesion and blood group antigen.

[Learn More](#)

Parameter	Value
Product Name	ERMAP Rabbit pAb
Remarks/Alias	RD; SC; BTN5; PRO2801; ERMAP
Species	Human
GeneID (Human)	114625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-155 of human ERMAP (NP_061008.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96PL5
Protein Weight	53kDa
Shipping	Ice bag