

Anti-Erp19 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O95881

Item Number: CM0006136



Introduction

High-quality rabbit monoclonal antibody against human ERp19 (thioredoxin domain-containing protein 12), validated for WB, IHC-P, and ELISA. Ideal for studying ER protein folding and disulfide bond formation.

[Learn More](#)

Parameter	Value
Product Name	ERp19 Rabbit mAb
Remarks/Alias	AG1; AGR1; ERP16; ERP18; ERP19; TLP19; hAG-1; PDIA16; hTLP19; ERp19
Species	Human
Gene ID (Human)	51060
Gene ID	51060
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ERp19 (NP_056997.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O95881
Protein Weight	19kDa
Shipping	Ice bag