

Anti-Erp29 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P30040

Item Number: CM0003621



Introduction

High-quality rabbit monoclonal antibody against human ERp29 (ERp28/ERp31), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ER protein folding and secretory pathway.

[Learn More](#)

Parameter	Value
Product Name	Erp29 Rabbit mAb
Remarks / Alias	ERp28; ERp31; PDIA9; PDI-DB; C12orf8; HEL-S-107; ERp29
Species	Human
GeneID (Human)	10961
GeneID	10961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human ERp29 (P30040).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30040
Protein Weight	29kDa
Shipping	Ice bag