

Anti-Erp29 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30040

Item Number: CM0028011



Introduction

Rabbit polyclonal antibody targeting human ERp29 (ERp28/ERp31). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ER protein folding and secretory pathway studies.

[Learn More](#)

Parameter	Specification
Product Name	ERp29 Rabbit pAb
Remarks/Alias	ERp28; ERp31; PDIA9; PDI-DB; C12orf8; HEL-S-107; ERp29
Species	Human
Gene ID (Human)	10961
Gene ID	10961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-261 of human ERp29 (NP_006808.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30040
Protein Weight	29kDa
Shipping	Ice bag