

Anti-Erp44 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bs26

Item Number: CM0026870



Introduction

ERP44 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with human, mouse, and rat. Useful for studying ER protein retention, oxidative folding, and calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	ERP44 Rabbit pAb
Remarks/Alias	PDIA10; TXNDC4; ERP44
Species	Human
Gene ID (Human)	23071
Gene ID	23071
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 57-406 of human ERP44 (NP_055866.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BS26
Protein Weight	47kDa
Shipping	Ice bag