

Anti-Erp57 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30101

Item Number: CM0014488



Introduction

Rabbit polyclonal antibody against ERp57 (PDIA3), validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. ERp57 is a protein disulfide isomerase and MHC I peptide-loading complex component. Ideal for studies on protein folding, redox regulation, and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	ERp57 Rabbit pAb
Remarks/Alias	P58; ER60; ERp57; ERp60; ERp61; GRP57; GRP58; PI-PLC; HsT17083; HEL-S-269; HEL-S-93n
Species	Human
GeneID (human)	2923
GeneID	2923
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 206-505 of human ERp57 (NP_005304.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30101
Protein Weight	54kDa
Shipping	Ice bag