

Anti-Erp57 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P30101

Item Number: CM0005635



Introduction

High-quality Anti-ERp57 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Targets protein disulfide-isomerase A3, a key chaperone in MHC I antigen presentation. Ideal for immunology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	ERp57 Rabbit mAb
Remarks/Alias	P58; ER60; ERp57; ERp60; ERp61; GRP57; GRP58; PI-PLC; Hst17083; HEL-S-269; HEL-S-93n
Species	Human
GeneID (human)	2923
GeneID	2923
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P30101
Protein Weight	57kDa
Shipping	Ice bag