

Anti-Ervw-1 (Syncytin-1) Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uqf0

Item Number: CM0020549



Introduction

Rabbit polyclonal antibody targeting human ERVW-1/Syncytin-1, validated for Western blot (WB) and ELISA. Detects the endogenous retroviral envelope protein involved in trophoblast fusion and multiple sclerosis pathogenesis.

[Learn More](#)

Parameter	Value
Product Name	ERVW-1 Rabbit pAb
Remarks/Alias	ENV; ENVW; HERVW; ERVWE1; HERV7Q; HERV-7q; HERVWENV; HERV-W-ENV; ERVW-1
Species	Human
GeneID (human)	30816
GeneID	30816
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 201-300 of human ERVW-1 (NP_055405.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UQF0
Protein Weight	60kDa
Shipping	Ice bag