

## Anti-Esam Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ap7

Item Number: CM0013825

### Introduction

Rabbit polyclonal antibody against human ESAM, validated for Western blot (WB) and ELISA. Ideal for endothelial cell adhesion, junctional biology, and vascular research. Predicted molecular weight: 41 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ESAM Rabbit pAb                                                                                                                       |
| Remarks/Alias    | W117m; ESAM                                                                                                                           |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 90952                                                                                                                                 |
| GeneID           | 90952                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-248 of human ESAM (NP_620411.2). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | Q96AP7                                                                                                                                |
| Protein Weight   | 41kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |