


## Anti-Esd (Esterase D) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P10768

Item Number: CM0028851



### Introduction

Rabbit polyclonal antibody against human Esterase D (S-formylglutathione hydrolase, ESD) for Western blot, IHC-P, and ELISA. Detects human and mouse ESD (31kDa) involved in formaldehyde detoxification.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ESD Rabbit pAb                                                                                                                      |
| Remarks/Alias    | FGH; ESD                                                                                                                            |
| Species          | Human                                                                                                                               |
| GeneID (Human)   | 2098                                                                                                                                |
| GeneID           | 2098                                                                                                                                |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human ESD (NP_001975.1). |
| Source           | Rabbit                                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                                    |
| Cross Reactivity | Human, Mouse                                                                                                                        |
| SWISS            | P10768                                                                                                                              |
| Protein Weight   | 31kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |