


Anti-Espl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14674

Item Number: CM0013476



Introduction

ESPL1 rabbit polyclonal antibody suitable for WB and ELISA applications. Recognizes human, mouse, and rat separase, a key mitotic protease that triggers chromosome separation by cleaving cohesin.

[Learn More](#)

Parameter	Value
Product Name	ESPL1 Rabbit pAb
Remarks/Alias	ESP1; SEPA; ESPL1
Species	Human
GeneID (Human)	9700
GeneID	9700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1220-1500 of human ESPL1 (NP_036423.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14674
Protein Weight	233kDa
Shipping	Ice bag