


Anti-Esrα Rabbit Polyclonal Antibody For Wb, Elisa - P03372

Item Number: CM0015078



Introduction

Rabbit polyclonal antibody targeting human estrogen receptor alpha (ESRα, ER-α). Suitable for Western blot (WB) and ELISA applications. Cross-reacts with human and mouse ESRα. For research use in hormone signaling, cancer biology, and drug development studies.

[Learn More](#)

Parameter	Value
Product Name	ESRα Rabbit pAb
Remarks/Alias	ER; ESR; Era; ESRA; ESTRR; NR3A1; ESRα
Species	Human
GeneID (Human)	2099
Immunogen	Recombinant fusion protein containing a sequence within amino acids 516-595 of human ESRα (NP_000116.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P03372
Protein Weight	66kDa
Shipping	Ice bag