


Anti-Esrra (Err-Alpha) Rabbit Polyclonal Antibody - P11474

Item Number: CM0013524



Introduction

Reliable ESRRR Rabbit Polyclonal Antibody for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat ESRRR. Supports energy metabolism and nuclear receptor studies.

[Learn More](#)

Parameter	Value
Product Name	ESRRR Rabbit pAb
Remarks/Alias	ERR1; ERRa; ESRL1; NR3B1; ERRalpha; ESRRR
Species	Human
GeneID (Human)	2101
GeneID	2101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human ESRRR (NP_004442.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11474
Protein Weight	46kDa
Shipping	Ice bag