


Anti-Esrra Rabbit Monoclonal Antibody For Wb, Elisa - P11474

Item Number: CM0005304



Introduction

High-quality rabbit monoclonal antibody targeting human ESRR (ERRalpha) for WB and ELISA applications, with cross-reactivity to human and mouse. Ideal for energy metabolism and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	ESRR Rabbit mAb
Remarks/Alias	ERR1; ERRa; ESRL1; NR3B1; ERRalpha; ESRR
Species	Human
GeneID (Human)	2101
GeneID	2101
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11474
Protein Weight	46kDa
Shipping	Ice bag