


Anti-Esrra Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11474

Item Number: CM0023990



Introduction

Rabbit polyclonal antibody against human ESRRRA (ERR-alpha), a nuclear receptor involved in energy metabolism. Validated for WB, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Suitable for transcriptional regulation and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	ESRRRA Rabbit pAb
Remarks/Alias	ERR1; ERRa; ESRL1; NR3B1; ERRalpha; ESRRRA
Species	Human
GeneID	2101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human ESRRRA (NP_004442.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11474
Protein Weight	46kDa
Shipping	Ice bag