


Anti-Esrrb Rabbit Polyclonal Antibody For Wb, Elisa - O95718

Item Number: CM0014741



Introduction

Rabbit polyclonal antibody against human ESRRB (ERR β) validated for Western blot and ELISA. Detects endogenous ESRRB in human, mouse, and rat samples. Ideal for stem cell pluripotency and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	ESRRB Rabbit pAb
Remarks/Alias	ERR2; ERRb; ESRL2; NR3B2; DFNB35; ERRbeta2; ERR beta-2; ESRRB
Species	Human
GeneID (Human)	2103
GeneID	2103
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ESRRB (NP_004443.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95718
Protein Weight	48kDa
Shipping	Ice bag