


Anti-Estrogen Receptor Alpha (Phospho S106) Rabbit Polyclonal Antibody For Wb, If/Icc - P03372

Item Number: CM0001855



Introduction

Phospho-specific rabbit polyclonal antibody targeting Estrogen Receptor alpha (ESR1) at serine 106. Reacts with human and mouse species. Suitable for WB and IF/ICC in cancer biology and endocrinology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ESRα-S106 Rabbit pAb
Remarks/Alias	ER; ESR; Era; ESRA; ESTRR; NR3A1; Phospho-ESRα-S106
Species	Human
GeneID (human)	2099
GeneID	2099
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Mouse
SWISS	P03372
Protein Weight	66kDa
Shipping	Ice bag