


Anti-Estrogen Receptor Beta Rabbit Polyclonal Antibody For Wb, Elisa - Q92731

Item Number: CM0030239



Introduction

Rabbit polyclonal antibody against human Estrogen Receptor beta (ER-beta, ESR2). Suitable for Western blot and ELISA, cross-reactive with mouse and rat. Targets nuclear hormone receptor involved in estrogen signaling in reproductive and other tissues. Ideal for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Estrogen Receptor beta Rabbit pAb
Remarks/Alias	Erb; ESRB; ODG8; ESTRB; NR3A2; ER-BETA; ESR-BETA; Estrogen Receptor beta
Species	Human
Gene ID	2100
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92731
Protein Weight	59kDa
Shipping	Ice bag