

Anti-Esyt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9BSj8

Item Number: CM0017666

Introduction

Rabbit polyclonal antibody against human ESYT1, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets extended synaptotagmin-1, an ER-cell membrane tether protein involved in calcium signaling and lipid transport. For research use.

[Learn More](#)

Parameter	Specification
Product Name	ESYT1 Rabbit pAb
Remarks/Alias	MBC2; FAM62A; ESYT1
Species	Human
Gene ID (Human)	23344
Gene ID	23344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-220 of human ESYT1 (NP_056107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
Swiss-Prot	Q9BSJ8
Protein Weight	123kDa
Shipping	Ice bag