

Anti-Etfa Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P13804

Item Number: CM0026840



Introduction

Rabbit polyclonal antibody against human electron transfer flavoprotein subunit alpha (ETFA), validated for WB, IHC-P, IP, and ELISA applications. Reacts with human, mouse, and rat. Ideal for studying mitochondrial fatty acid oxidation and glutaricaciduria type II research.

[Learn More](#)

Parameter	Value
Product Name	ETFA Rabbit pAb
Remarks/Alias	EMA; GA2; MADD; ETFA
Species	Human
Gene ID (Human)	2108
Gene ID	2108
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-333 of human ETFA (NP_000117.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13804
Protein Weight	35kDa
Shipping	Ice bag