

Anti-Etfa Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P13804

Item Number: CM0010918



Introduction

Rabbit monoclonal antibody against ETFA, suitable for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets mitochondrial electron transfer flavoprotein, essential for fatty acid oxidation and linked to glutaric acidemia type II.

[Learn More](#)

Parameter	Value
Product Name	ETFa Rabbit mAb
Remarks/Alias	EMA; GA2; MADD
Species	Human
GeneID (Human)	2108
GeneID	2108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-333 of human ETFa (NP_000117.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13804
Protein Weight	35kDa
Shipping	Ice bag