

Anti-Etfa Rabbit Polyclonal Antibody [Ko Validated] - P13804

Item Number: CM0023696



Introduction

Rabbit polyclonal antibody against ETFA, validated by knockout. Suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, rat ETFA. Ideal for mitochondrial fatty acid oxidation and amino acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ETFA Rabbit pAb
Remarks/Alias	EMA; GA2; MADD; FA
Species	Human
GeneID (human)	2108
GeneID	2108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-333 of human ETFA (NP_000117.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13804
Protein Weight	35kDa
Shipping	Ice bag