

Anti-Etfdh Monoclonal Antibody For Wb, If, Ihc, Elisa - Q16134

Item Number: CM0011513

Introduction

Validated rabbit monoclonal antibody against human ETFDH (MADD, ETFQO), suitable for WB, IF/ICC, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat ETFDH protein, a key enzyme in mitochondrial electron transport.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ETFDH Rabbit mAb
Remarks/Alias	MADD; ETFQO
Species	Human
GeneID	2110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-617 of human ETFDH (NP_004444.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16134
Protein Weight	68kDa
Shipping	Ice bag