

Anti-Etfdh Rabbit Polyclonal Antibody For Wb And Elisa - Q16134

Item Number: CM0030515

Introduction

Rabbit polyclonal antibody targeting human ETFDH (electron transfer flavoprotein-ubiquinone oxidoreductase), validated for Western blot and ELISA. Detects human, mouse, rat. Research use in fatty acid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ETFDH Rabbit pAb
Remarks/Alias	MADD; ETFQO; ETFDH
Species	Human
Gene ID (Human)	2110
Gene ID	2110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 318-617 of human ETFDH (NP_004444.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16134
Protein Weight	68kDa
Shipping	Ice bag