


Anti-Ethe1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95571

Item Number: CM0013427



Introduction

ETHE1 Rabbit pAb detects human, mouse, rat ETHE1. Suitable for WB, IHC-P, IF/ICC, ELISA. ETHE1 is a mitochondrial persulfide dioxygenase essential for hydrogen sulfide catabolism.

[Learn More](#)

Parameter	Value
Product Name	ETHE1 Rabbit pAb
Remarks/Alias	HSCO; YF13H12; ETHE1
Species	Human
GeneID (Human)	23474
GeneID	23474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 8-245 of human ETHE1 (NP_055112.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95571
Protein Weight	28kDa
Shipping	Ice bag