


Anti-Ets1 Polyclonal Antibody For Wb, If-P, Elisa - P14921

Item Number: CM0015130



Introduction

Rabbit polyclonal antibody against human ETS1 transcription factor, validated for WB, IF-P, ELISA, and reactive with human, mouse, rat. Ideal for studying lymphoid differentiation, angiogenesis, and cytokine regulation.

[Learn More](#)

Parameter	Specification
Product Name	ETS1 Rabbit pAb
Remarks/Alias	p54; ETS-1; EWSR2; c-ets-1; ETS1
Species	Human
GeneID (human)	2113
GeneID	2113
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 342-441 of human ETS1 (NP_005229.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14921
Protein Weight	50kDa
Shipping	Ice bag