


Anti-Ets1 Rabbit Monoclonal Antibody For Wb And Elisa - P14921

Item Number: CM0003657



Introduction

High-quality rabbit monoclonal antibody targeting human ETS1 (p54), validated for Western blot and ELISA. Suitable for transcription factor studies, lymphoid cell differentiation, and angiogenesis research. Reference: P14921.

[Learn More](#)

Parameter	Value
Product Name	ETS1 Rabbit mAb
Remarks/Alias	p54; ETS-1; EWSR2; c-ets-1
Species	Human
GeneID (human)	2113
GeneID	2113
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ETS1 (P14921).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14921
Protein Weight	50kDa
Shipping	Ice bag