

## Anti-Ets2 Polyclonal Antibody For Wb, If/Icc, Elisa - P15036

Item Number: CM0030212



### Introduction

Rabbit polyclonal antibody against human ETS2, suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for transcription regulation and cancer research studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ETS2 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | ETS2IT1; ETS2                                                                                                                         |
| Species          | Human                                                                                                                                 |
| GeneID (human)   | 2114                                                                                                                                  |
| GeneID           | 2114                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 76-170 of human ETS2 (NP_005230.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | P15036                                                                                                                                |
| Protein Weight   | 53kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |