


## Anti-Etv1 Monoclonal Antibody For Wb, Ip, Elisa, Chip - P50549

Item Number: CM0005177



### Introduction

Rabbit monoclonal antibody targeting human ETV1 (ER81), validated for WB, IP, ELISA, and ChIP. Cross-reacts with mouse and rat. Suitable for transcription factor and oncology studies.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] ETV1 Rabbit mAb                                                                      |
| Remarks/Alias    | ER81                                                                                                |
| Species          | Human                                                                                               |
| GeneID (Human)   | 2115                                                                                                |
| GeneID           | 2115                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ETV1 (NP_004947.2). |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibodies                                                                               |
| Application      | WB, IP, ELISA, ChIP                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                   |
| SWISS            | P50549                                                                                              |
| Protein Weight   | 55 kDa                                                                                              |
| Shipping         | Ice bag                                                                                             |