


Anti-Etv1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P50549

Item Number: CM0013512



Introduction

ETV1 rabbit polyclonal antibody validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying ETV1 transcriptional regulation, neurodifferentiation, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	ETV1 Rabbit pAb
Remarks/Alias	ER81; ETV1
Species	Human
Gene ID (Human)	2115
Gene ID	2115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-300 of human ETV1 (NP_004947.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50549
Protein Weight	55kDa
Shipping	Ice bag