

Anti-Etv1/Er81 Rabbit Polyclonal Antibody For Wb & Elisa - P50549

Item Number: CM0025109



Introduction

Rabbit polyclonal antibody specific for human ETV1/ER81. Validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for studying transcription, neuronal differentiation, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	ETV1 Rabbit pAb
Remarks/Alias	ER81; ETV1
Species	Human
Gene ID (Human)	2115
Gene ID	2115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ETV1 (NP_004947.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50549
Protein Weight	55kDa
Shipping	Ice bag