


Anti-Etv4 Polyclonal Antibody For Wb And Elisa - P43268

Item Number: CM0013898



Introduction

CamelBio rabbit polyclonal antibody against human ETV4 (E1AF/PEA3) for WB and ELISA, cross-reacts with mouse and rat. Ideal for transcription and neurogenesis studies, supporting motor neuron and keratinocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	ETV4 Rabbit pAb
Remarks/Alias	E1AF; PEA3; E1A-F; PEAS3; ETV4
Species	Human
GeneID (human)	2118
GeneID	2118
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 79-178 of human ETV4 (NP_001977.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43268
Protein Weight	54kDa
Shipping	Ice bag