

Anti-Etv5 Polyclonal Antibody For Wb And Elisa - P41161

Item Number: CM0022628



Introduction

Anti-ETV5 rabbit polyclonal antibody for WB and ELISA applications. Recognizes human, mouse, rat ETV5/ERM protein (58 kDa), a DNA-binding transcription factor involved in gene regulation. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	ETV5 Rabbit pAb
Remarks/Alias	ERM; ETV5
Species	Human
GeneID (Human)	2119
GeneID	2119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-510 of human ETV5 (P41161).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41161
Protein Weight	58kDa
Shipping	Ice bag