

Anti-Etv6 Polyclonal Antibody For Wb, Ihc-P, Elisa - P41212

Item Number: CM0022680



Introduction

Rabbit polyclonal antibody against human ETV6 (TEL) for WB, IHC-P, ELISA. Validated in human, mouse, rat. ETV6 is a transcriptional repressor involved in hematopoiesis and malignant transformation. Ideal for cancer and transcription regulation studies.

[Learn More](#)

Parameter	Description
Product Name	ETV6 Rabbit pAb
Remarks/Alias	TEL; THC5; TEL/ABL; ETV6
Species	Human
GeneID (Human)	2120
GeneID	2120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ETV6 (NP_001978.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41212
Protein Weight	53kDa
Shipping	Ice bag