

Anti-Etv6 Rabbit Polyclonal Antibody For Wb, Elisa - P41212

Item Number: CM0016566



Introduction

High-quality rabbit polyclonal antibody targeting human ETV6 (TEL), a transcriptional repressor validated for WB and ELISA. Cross-reacts with mouse and rat; 53kDa. Ideal for hematopoiesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ETV6 Rabbit pAb
Remarks/Alias	TEL, THCS, TEL/ABL, ETV6
Species	Human
Gene ID (Human)	2120
Gene ID	2120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ETV6 (NP_001978.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41212
Protein Weight	53kDa
Shipping	Ice bag