


Anti-Etv7 Polyclonal Antibody For Wb, Elisa - Q9Y603

Item Number: CM0026032



Introduction

ETV7 Rabbit pAb is a polyclonal antibody for human ETV7 detection in WB, ELISA. Validated for research use, targeting the transcription factor ETV7 involved in hematopoietic gene regulation.

[Learn More](#)

Parameter	Value
Product Name	ETV7 Rabbit pAb
Remarks/Alias	TEL2; TELB; TEL-2; ETV7
Species	Human
GeneID (Human)	51513
GeneID	51513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human ETV7 (NP_057219.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y603
Protein Weight	39kDa
Shipping	Ice bag