

Anti-Eva1A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H8M9

Item Number: CM0027882

Introduction

Rabbit polyclonal antibody against human EVA1A (FAM176A/TMEM166) for WB, IHC-P, ELISA. Regulates autophagy and apoptosis. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	EVA1A Rabbit pAb
Remarks/Alias	FAM176A; TMEM166; EVA1A
Species	Human
GeneID (Human)	84141
GeneID	84141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-152 of human EVA1A (NP_115557.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8M9
Protein Weight	17kDa
Shipping	Ice bag